

AestheLine

Sterile single use polydioxanone suture with needle





AestheLine MONO

- ✓ The most common type of thread
- ✓ Made from PDO (also for some sizes can be made from PLLA, PCL)
- ✓ Non cog
- ✓ Easy to set up

INDICATIONS

- Used for loss of skin tone of the face and body
- Additional enhancement of the effectiveness of the procedure when setting other types of threads
- Maintaining the effect after the thread lifting procedure
- Reduced tone and presence of wrinkles in areas with thin tissues
 - The area of the outer corner of the eye, around the mouth, neck area
- Reduced tone and the presence of wrinkles in small areas
 - The central line of the neck, the area between the eyebrows, the outer corner of the eye

THREAD DETAILS




PACKAGE - 20ea (10ea x 2 EO pouches)

CODE	N/G	N/L	T/L	USP	PDO	PLLA	PCL
MONO 31G 12mm	31G	12	16	7/0	●	●	
MONO 31G 25mm		25	30		●	●	●
MONO 30G 13mm	30G	13	15	6/0	●	●	●
MONO 30G 19mm		19	16		●	●	●
MONO 30G 25mm		25	30		● CE	●	●
MONO 30G 38mm		38	50		● CE	●	●
MONO 29G 25mm	29G	25	30	6/0	● CE	●	●
MONO 29G 38mm		38	50		● CE	●	●
MONO 29G 50mm		50	70		● CE	●	●
MONO 29G 60mm		60	80		● CE	●	●
MONO 27G 38mm	27G	38	50	5/0	● CE	●	●
MONO 27G 50mm		50	70		● CE	●	●
MONO 27G 60mm		60	80		● CE	●	●
MONO 27G 60mm		60	90		● CE	●	●
MONO 27G 90mm	26G	90	140	5/0	● CE	●	●
MONO 26G 38mm		38	50		● CE	●	
MONO 26G 50mm		50	70		● CE	●	
MONO 26G 60mm		60	80		● CE	●	
MONO 26G 70mm		70	90		●	●	
MONO 26G 90mm		90	150		● CE	●	

AestheLine SCREW

- ✓ Spiral yarn configuration
- ✓ It looks like a spiral, made from PDO (also for some sizes can be made from PLLA, PCL)
- ✓ Provides a pronounced tissue alteration, a prolonged physiological inflammatory response and, as a consequence, a pronounced dermorestructuring effect
- ✓ Indications for use are the same as for monofilaments, but due to the spiral structure they have longer biodegradation times, provide a more pronounced tissue compaction and a prolonged effect

THREAD DETAILS



PACKAGE - 20ea (10ea x 2 EO pouches)

CODE	N/G	N/L	T/L	USP	PDO	PLLA	PCL
SCREW 31G 25mm	31G	25	30	7/0	●	●	●
SCREW 30G 25mm	30G	25	30	6/0	● CE	●	●
SCREW 30G 38mm		38	50		● CE	●	●
SCREW 29G 25mm	29G	25	30	6/0	● CE	●	●
SCREW 29G 38mm		38	50		● CE	●	●
SCREW 29G 50mm		50	70		● CE	●	●
SCREW 29G 60mm		60	80		● CE	●	●
SCREW 27G 50mm	27G	50	70	5/0	● CE	●	●
SCREW 27G 60mm		60	80		● CE	●	●
SCREW 26G 60mm	26G	60	80	5/0	● CE	●	●

AestheLine MULTI

- ✓ Flat single threads with a loose free end in the form of a "tassel"
- ✓ Non cog, PGA (Polyglycolic Acid) material
- ✓ Each fine thread, united into a common structure, after installation in the tissue contributes to a more pronounced restructuring of the skin

INDICATIONS

- Stretch marks in various areas of the body after pregnancy, fluctuations in its weight
- Scars after acne
- Post-surgical scars
- Post-traumatic scars
- Furrows and creases on the skin of the face and body - "Rings of Venus" (horizontal creases on the neck), skin creases in the projection of the nasolabial fold

THREAD DETAILS



PACKAGE - 20ea (10ea x 2 EO pouches)

CODE	N/G	N/L	T/L	USP	PGA
MULTI PGA 29G 25mm	29G	25	30	6/0	●
MULTI PGA 29G 38mm		38	50		●
MULTI PGA 29G 50mm		50	70		●
MULTI PGA 29G 60mm		60	80		●
MULTI PGA 27G 38mm	27G	38	50	5/0	●
MULTI PGA 27G 50mm		50	70		●
MULTI PGA 27G 60mm		60	80		●

AestheLine TWIN

- ✓ 2 threads are twisted
- ✓ Non cog, made of PDO
- ✓ Create effective volume

THREAD DETAILS



PACKAGE - 20ea (10ea x 2 EO pouches)

CODE	N/G	N/L	T/L	USP	PDO
TWIN 26G 38mm	26G	38	50	7/0	●
TWIN 26G 60mm		60	80		●
TWIN 26G 90mm		90	140		●

AestheLine


TWIN SCREW

- ✓ Two threads are twisted together and twisted into a spiral around the needle
- ✓ Non cog, made of PDO
- ✓ After the needle with thread is set in the fabric, the guide needle is pulled out, and the two twisted threads begin to unravel, and thus additional volume is given. Wrinkles are instantly smoothed out due to the formed frame of the established threads

INDICATIONS

- For the correction of wrinkles in areas with mobile facial expressions
- For severe expression lines (puppets, nasolabial folds, wrinkles in the mouth and eyes)
- A spring thread introduced into a wrinkle instantly smoothes it out, while keeping facial expressions natural

THREAD DETAILS





PACKAGE - 20ea (10ea x 2 EO pouches)

CODE	N/G	N/L	T/L	USP	PDO
TWIN SCREW 26G 38mm	26G	38	50	7/0	●
TWIN SCREW 26G 50mm		50	70		●
TWIN SCREW 26G 60mm		60	80		●
TWIN SCREW 26G 90mm		90	140		●

AestheLine BOOM-6

- ✓ Multicomponent thread. Consisting of 6 fine structured threads
- ✓ Non cog, made by PDO (could be made by PLLA also)
- ✓ Creates effective volume (filling nasolabial folds and wrinkles)

THREAD DETAILS



PACKAGE - 20ea (4ea x 5 EO pouches)

CODE	N/G	N/L	T/L	USP	PDO	PLLA
BOOM-6 23G 38mm	23G	38	50	6/0 6ea	● CE	●
BOOM-6 23G 60mm		60	80		● CE	●

AestheLine BOOM-7

- ✓ Multicomponent thread. Consisting of 7 fine structured threads
- ✓ Non cog, made by PDO (could be made PLLA also)
- ✓ Creates effective volume (filling nasolabial folds and wrinkles)

THREAD DETAILS



PACKAGE - 20ea (4ea x 5 EO pouches)

CODE	N/G	N/L	T/L	USP	PDO
BOOM-7 23G 38mm	23G	38	50	5/0 7ea	●
BOOM-7 23G 60mm		60	80		● CE

AestheLine BOOM-10

- ✓ Multicomponent thread. Consisting of 10 fine structural threads
- ✓ Non cog, made by PDO (could be made of PLLA also)
- ✓ Create extremely effective volume

INDICATIONS

- Elimination of folds on the face: nasolabial, "puppet" lines if you do not want to resort to filler injections (for example, if you have a tendency to swelling) or to improve the results after their application
- Thickening of the skin on the face and body in areas of maximum laxity and reduced tone

THREAD DETAILS



PACKAGE - 20ea (4ea x 5 EO pouches)

CODE	N/G	N/L	T/L	USP	PDO
BOOM-10 23G 38mm	23G	38	50	7/0 10ea	●

AestheLine

N

- ✓ Thin but it's exquisite cog thread
- ✓ Use for eye area
- ✓ Made of PDO
- ✓ The outer part of the thread, twisted in a spiral, has multidirectional notches, which provides a pronounced stimulating effect on the tissue. The inner part of the thread, installed in the guide (thin, sharp needle), has converging notches and is divided into two halves : lifting and fixing, which provides maximum effect on tissues

INDICATIONS

- Loss of volume in the zygomatic region with pre-existing ptosis and decreased tone
- Loss of tone in areas with dynamic activity - Buccal region

THREAD DETAILS



PACKAGE - 20ea (4ea x 5 EO pouches)

CODE	N/G	N/L	T/L	USP	PDO
N 25G 38mm	25G	38	136	5/0	●
N 25G 60mm		60	185		●

AestheLine CLL

- ✓ Long thread inside and outside the conductor notches over the entire plane, 360 degrees, due to which it holds and fixes tissues for a long time
- ✓ With cannula L type
- ✓ Good for making face line

THREAD DETAILS



PACKAGE - 20ea (2ea x 10 EO pouches)

CODE	N/G	N/L	T/L	USP	PDO
CLL 18G 100mm	18G	100	230	2	●
CLL 19G 100mm	19G	100	230	0	●
CLL 21G 60mm	21G	60	150	2/0	●
CLL 21G 90mm		90	210		●
CLL 23G 60mm	23G	60	150	3/0	●
CLL 23G 90mm		90	210		●

AestheLine M FILL

- ✓ Can be used with filler
- ✓ Effective for nasolabial folds
- ✓ Made of PDO
- ✓ The thread is wound in a tight spiral on a safe atraumatic guide. A special thread insertion aid helps it to easily overcome tissue resistance without losing its spiral structure. Due to its special configuration, the thread stimulates tissue several times more intensively and induces collagen formation

INDICATIONS

- Dropped corners of the mouth
- Lack of volume in the zygomatic zone if you do not want to resort to filler injections (for example, if you have a tendency to swelling, the presence of hernias) or to improve their results
- Loss of definition of facial contours
- Pronounced furrows on the face - for example, cheekbone

THREAD DETAILS



PACKAGE - 20ea (2ea x 10 EO pouches)

CODE	N/G	N/L	T/L	USP	PDO
M FILL 21G 60mm	21G	60	310	3/0	●
M FILL 21G 60mm			350		●
M FILL 21G 60mm			410		●
M FILL 21G 60mm			450	4/0	●
M FILL 21G 60mm			530		●
M FILL 21G 60mm			600	5/0	●
M FILL 23G 60mm	23G	60	360	4/0	●
M FILL 23G 60mm			410		●
M FILL 23G 60mm			450	5/0	●
M FILL 23G 60mm			530		●
M FILL 25G 50mm	25G	50	320	5/0	●
M FILL 25G 50mm			380		●
M FILL 27G 50mm	27G	50	270	5/0	●
M FILL 29G 50mm	29G	50	290	6/0	●
TM FILL 21G 60mm	21G	60	150	4/0	●

AestheLine

COG MESH


- ✓ Combination yarn, which combines mesh and cog.
(A thread with a notch is inserted inside the mesh thread)
- ✓ The mesh thread fills in wrinkles, and at the same time, the knurled thread tightens the fabrics
- ✓ Made of PDO and L cannula or W cannula
- ✓ The thread insertion procedure is the same as for the cannula-type serrated thread

THREAD DETAILS


CW




CL



CW



CL



PACKAGE - 20ea (2ea x 10 EO pouches)

CODE	N/G	N/L	T/L	USP	PDO
COG MESH CW 18G 100mm	18G	100	155	2	●
COG MESH CL 18G 100mm		100	155		●


AestheLine MESH COVER

- ✓ Made with mesh type filaments
- ✓ Non cog, made of PCL
- ✓ PDO thread, consisting of 8 woven PCL threads 6.0 in the form of a hollow tube
- ✓ Strong volumizing ability
- ✓ Excellent drainage, which is used in the correction of nasolabial and labial folds, restoration of the oval line along the edge of the lower jaw and a sharper mandibular line
- ✓ The ability to use in conjunction with filler

INDICATIONS

- Deep nasolabial folds
- Deep lip-chin folds
- Restoration of the oval line along the edge of the lower jaw
- Creating a clearer mandibular line
- Reduction of wings

THREAD DETAILS



PACKAGE - 20ea (2ea x 10 EO pouches)

CODE	N/G	N/L	T/L	USP	PCL
MESH COVER PCL CW 21G 60mm	21G	60	25	6/0	●
MESH COVER PCL CW 21G 60mm			30		●
MESH COVER PCL CW 21G 60mm			40		●
MESH COVER PCL CW 23G 60mm	23G	60	25	6/0	●
MESH COVER PCL CW 23G 60mm			30		●
MESH COVER PCL CW 23G 60mm			40		●


AestheLine RHINO

- ✓ Specialized thread for elevation of nose tip and bridge
- ✓ Made of PDO and also could be PLLA (PCL also for some size)

INDICATIONS

- Deformation of the dorsum of the nose ("saddle" deformity, etc.)
- Patient's desire to visually narrow the nose
- Dropped nose tip

THREAD DETAILS



PACKAGE - 20ea (2ea x 10 EO pouches)

CODE	N/G	N/L	T/L	USP	PDO	PLLA	PCL
R 19G 38mm	19G	38	40	2	● CE	●	●
R 19G 38mm			50		● CE	●	●
R 19G 50mm		50	60		● CE	●	●
R 19G 60mm		60	80		● CE	●	●
R 21G 60mm	21G	60	80	2/0	● CE	●	●

AestheLine






DOUBLE ARM

- ✓ Strong, bulky cast thread with multiple three-dimensional fixation mechanisms in rocket-shaped fabrics on two carrier needles. The innovative configuration of the thread creates a large area of contact of its surface with the tissues, pronounced fixation and enhanced fibrotic response, provides the most pronounced and long-lasting lifting of facial tissues, comparable to surgical methods of aesthetic correction
- ✓ Both side of thread, there are attached 19G or 21G needle (No need hook for procedure)
- ✓ Our best recommendation thread
- ✓ Made of PDO (PCLLA also)

INDICATIONS

- Loss of volume and pronounced age-related changes in the middle and lower third of the face
- Significant loss of definition of facial contours, presence of "heavy flaps" and double chin
- The need to obtain the most pronounced and lasting correction of problem areas

THREAD DETAILS



PACKAGE - 2ea (1ea x 2 EO pouches)

CODE	N/G	N/L	T/L	USP	PDO	PLLA	PCL	PCLLA
LN 19G 50mm	19G	50	200	0	•	•	•	•
LN 19G 100mm		100	300		•	•	•	•
LN 19G 120mm		120	420		•	•	•	•
LN 21G 100mm	21G	100	300	0	•	•	•	•
LN 21G 120mm		120	400		•		•	•


AestheLine LN M1

- ✓ Same indications as Double arm, but with pressed cone-shaped incisions for movement and prolonged fixation of heavy integumentary tissues

INDICATIONS

- Dense heavy integumentary tissues of the face

THREAD DETAILS



PACKAGE - 2ea (1ea x 2 EO pouches)

CODE	N/G	N/L	T/L	USP	PDO
LN M1 19G 100mm	19G	100	320	0	●
LN M1 19G 100mm				1	●
LN M1 19G 120mm		120	320	0	●
LN M1 19G 120mm				1	●

AestheLine

EYE CANNULA

- ✓ Non cog, thread with U shape (edge) cannula
- ✓ Painless injection
- ✓ Made of PDO (could be made of PLLA, PCL also)
- ✓ A thread designed specifically for the eye area and located in an atraumatic guide in the form of a cannula. Allows you to correct the most delicate area on the face around the eyes without hematomas and rehabilitation period

INDICATIONS

- "Crow's feet" - wrinkles of the outer corner of the eye
- Painful hernias - increased fatty accumulations in the infraorbital region
- Decreased skin tone around the eyes

THREAD DETAILS



PDO



PLLA



PCL

PACKAGE - 20ea (4ea x 5 EO pouches)

CODE	N/G	N/L	T/L	USP	PDO	PLLA	PCL
E 30G 25mm	30G	25	35	6/0	● CE	●	●
E 30G 38mm		38	50		● CE	●	●


AestheLine CANNULA L

- ✓ Cog thread with L shape(edge) cannula
- ✓ Painless injection
- ✓ Made of PDO (could be made of PLLA, PCL also)
- ✓ The thread has 360-degree notches arranged in a spiral, thanks to which it holds and fixes tissues for a longer time

INDICATIONS

- Pronounced signs of gravitational ptosis - an influx of tissues above the nasolabial fold, loss of clarity of facial contours, "double chin", presence of "flew"
- Dense heavy integumentary tissues of the face

THREAD DETAILS



PACKAGE - 20ea (2ea x 10 EO pouches)

CODE	N/G	N/L	T/L	USP	PDO	PLLA	PCL
CL 18G 100mm	18G	100	155	2	● CE	●	●
CL 19G 100mm	19G	100	155	0	● CE	●	
CL 21G 60mm	21G	60	115	2/0	● CE	●	
CL 21G 90mm		90	165		● CE	●	
CL 23G 60mm	23G	60	115	3/0	● CE	●	
CL 23G 90mm		90	165		● CE		

AestheLine CANNULA W

- ✓ The thread has 360-degree notches arranged in a spiral, thanks to which it holds and fixes tissues for a longer time

INDICATIONS

- Pronounced signs of gravitational ptosis - an influx of tissues above the nasolabial fold, loss of clarity of facial contours, "double chin", presence of "flew"
- Dense heavy integumentary tissues of the face

THREAD DETAILS



PDO



PLLA



PCL

PACKAGE - 20ea (2ea x 10 EO pouches)

CODE	N/G	N/L	T/L	USP	PDO	PLLA	PCL
CW 18G 100mm	18G	100	150	2	● CE	●	●
CW 19G 100mm	19G	100	150	1	● CE	●	●
CW 21G 60mm	21G	60	135	2/0	● CE	●	●
CW 21G 90mm		90	145		● CE	●	●

AestheLine



M1

- ✓ Cast Thread with pressed conical notches for movement and prolonged fixation of heavy integumentary tissues

INDICATIONS

- Loss of volume and pronounced tissue prolapse in the zygomatic zone
- Significant loss of definition of facial contours
- The presence of "heavy wings"
- Pronounced volume of the "double chin"
- Severe loss of clarity of the cervico-chin angle
- Desire to obtain the longest possible correction result

THREAD DETAILS




PACKAGE - 20ea (2ea x 10 EO pouches)

CODE	N/G	N/L	T/L	USP	PDO	PCLLA
M1 CW 18G 100mm	18G	100	170	3&4	• CE	•
M1 CL 18G 100mm			175		• CE	•

AestheLine M2

- ✓ Thread with pressed notches, cannula type W
- ✓ Made according to the same principle as M1, but has a different shape of notches. This shape holds heavy skin better
- ✓ Effectively tightens the entire face and contour
- ✓ The insertion procedure is the same as for other cannula-type serrated suture
- ✓ Made of PDO

THREAD DETAILS



PACKAGE - 20ea (2ea x 10 EO pouches)

CODE	N/G	N/L	T/L	USP	PDO
M2 CW 18G 100mm	18G	100	170	3&4	● CE
M2 CW 19G 100mm	19G	100	170		● CE

AestheLine G LINE

- ✓ Made of 99.99% Gold
- ✓ Guide - sharp needle
- ✓ Effectively removes wrinkles around the eyes and whitens the face
- ✓ The procedure is the same as Mono or Screw

THREAD DETAILS



PACKAGE - 20ea (10ea x 2 EO pouches)

CODE	N/G	N/L	T/L	USP	GOLD
G LINE 29G 25mm	29G	25	32.5	5/0	●
G LINE 29G 38mm		38	50		●
G LINE 30G 13mm	30G	13	15	5/0	●
G LINE 30G 25mm		25	32.5		●
G LINE 30G 38mm		38	50		●

AestheLine TRIPLE

- ✓ Three threads are twisted together, on a sharp needle
- ✓ Non cog, made of PDO
- ✓ Create extremely effective volume

THREAD DETAILS



PACKAGE - 20ea (4ea x 5 EO pouches)

CODE	N/G	N/L	T/L	USP	PDO
TRIPLE 26G 38mm	26G	38	50	7/0	●
TRIPLE 26G 50mm		50	70		●
TRIPLE 26G 60mm		60	80		●
TRIPLE 26G 90mm		90	140		●

AestheLine

BLUNT NEEDLE

- ✓ Used blunt needle and cog thread
- ✓ Painless injection and less bleeding and swelling than sharp needle
- ✓ Made of PDO

THREAD DETAILS



PACKAGE - 20ea (4ea x 5 EO pouches)

CODE	N/G	N/L	T/L	USP	PDO
BLUNT NEEDLE 19G 100mm	19G	100	150	0	● CE

AestheLine CSL

- ✓ It is made by PDO with both side cog (inside and outside of needle)
- ✓ With sharp needle
- ✓ It is good for Making face line

INDICATIONS

- Use 1 thread, and get two threads effect
- Normal cogs use USP 5/0 and they are weak, but on CSL USP 3/0 is used, thus cog power is very strong
- It is multi-directional cog, so it has good lifting effect
- Can be used on overall face line, also thin skin layer such as eye line and forehead. But as it is sharp needle, not cannula, must be very careful while treating eye area or thin layers

THREAD DETAILS



PACKAGE - 20ea (4ea x 5 EO pouches)

CODE	N/G	N/L	T/L	USP	PDO
CSL 23G 38mm	23G	38	188	3/0	●
CSL 23G 60mm		60	240		●

Competitive
price

Long
experience

Special
skill