





Mesotherapy needles are very thin - this is done in order to reduce pain during the procedure. Outside diameter size is indicated by the letter G (on the Gage scale). As a rule, meso needles with sizes from 25 to 34G are used for meso, the most popular are 30 and 32G. This allows the

They consist of a steel tube, one of the ends is cut at an angle, and the other is fixed with a

special glue in a plastic base.

mesotherapeutic agent to be introduced to a shallow depth and excludes the occurrence of The length of the meso needles depends on the area to be treated and ranges from 3 to 25 mm. Mesotherapy of the face involves the use of short needles, and the use of long needles for the body. In addition, mesotherapy needles with corrective length are also produced.



N/G

N/L

13mm

13mm

4mm

Quantity

HYPODERMIC NEEDLE 25G 13mm 25G **HYPODERMIC NEEDLE** 27G 13mm 27G

CODE

HYPODERMIC NEEDLE 30G 4mm

👸 BienMESO

HYPODERMIC NEEDLE 30G 6mm	30G	6mm	
HYPODERMIC NEEDLE 30G 13mm		13mm	
HYPODERMIC NEEDLE 31G 4mm		4mm	
HYPODERMIC NEEDLE 31G 6mm	31G	6mm	10ea
HYPODERMIC NEEDLE 31G 13mm		13mm	
HYPODERMIC NEEDLE 32G 4mm	32G	4mm	
HYPODERMIC NEEDLE 32G 6mm		6mm	
HYPODERMIC NEEDLE 33G 4mm	33G	4mm	
HYPODERMIC NEEDLE 34G 4mm	34G	4mm	
⊯ Bien€	^ANN	JIII Δ	



The nanocannulas are made of stainless steel (SUS34) using micronized ultra-fine grinding technology, making them flexible and durable.



 Ultra-thin walls and micro-polishing of their surface - less tissue resistance during the procedure Widened inner diameter for insertion of the tightest fillers Marking along the length of the cannula to control the depth of its insertion into the tissue

Benefits of BienCannula

BienCannula is used for

Facial contouring

Bio-reinforcement Correction of wrinkles Creation of stiffeners Sculptural face correction

Volumetric plastic (volumetry)

BienCannula Type

CODE

GENERAL PURPOSE CANNULA 14G 100mm

GENERAL PURPOSE CANNULA 17G 100mm	17G	100m
GENERAL PURPOSE CANNULA 22G 50mm	22G	50mr
GENERAL PURPOSE CANNULA 22G 70mm	226	70mr
GENERAL PURPOSE CANNULA 25G 50mm	25G	50mr



N/G

14G

N/L

100mm

Quantity

20ea

